

# Agilent 7890 GC and 7000B GC/MS/MS

## Presentation for CA Dept of Toxic Substances Control

### Environmental Chemistry Laboratory

Jan 26, 2011



Agilent Technologies  
Tim Borrego  
Account Manager  
Santa Clara, CA



# Agenda

- Why MS/MS?
- Hardware Features of 7000B MS/MS
- Hardware Features of 7890 GC that help mass spec
- Software
- Applications



# Why MS/MS?

## Greater Selectivity Than SIM

### EI-SIM

Selectivity proportional to spectral resolution

No selectivity against ions with same  $m/z$


*unit mass resolution*

interference

analyte

### EI-MS/MS

Precursor selectivity same as SIM  
High probability that at least one product ion will be a unique dissociation product of the precursor **BUT** not the interference




The precursor ion should **NOT** be used for ion ratios or quantitation since the interferences will be the same as the SIM ion


# Why MS/MS?

## Lower detection limits by reducing noise


El: many ions from the source



Q1 SIM isolate precursor before CID





Product ion signal often decreases, but the percentage decrease in noise is much larger for real samples; S/N and detection limits improve



# Why GC/MS/MS?

## A Picture Is Worth a Thousand Words



GC/MS Single Quad SIM

GC/MS/MS MRM

EI 100fg HCB in “DIRTY”  
Matrix

*A chromatographer's  
dream: single peak on  
flat baseline*

Agilent took the best technologies from their industry leading 5975 GC/MSD:

- Heated monolithic gold plated quartz quadrupole
- Proven reliable high performance source design
- AUTOTUNE

and the 6400 Series LC/MS/MS:

- Linear acceleration enhanced Collision Cell
- Wide Mass-Bandwidth MS/MS ion optics
- Mass Hunter software



# Why Heated Quartz “Gold” Quads?


Unlike LC/MS, many high boiling neutral molecules enter the source and manifold of a GC/MS or GC/MS/MS

- Higher temperature reduces potential for contamination
  - 200C max
  - Virtually eliminates the need to clean quads
- Low coefficient for thermal expansion for quartz
  - Stable structure during maintenance cycles (hot-cool-hot) for source or detector

More stable tunes and methods over a longer period of time in real world sample environments




# GC/MS/MS Analyzer






# GC/MS/MS Analyzer



# Analyzer Doors Open



# Proven Inert Source Performance

Stay-Clean design extends maintenance intervals with dirty samples

Can be heated to 350°C

Has both repeller and extraction technologies to increase ion flux into the mass analyzer

Source tune parameters in the Autotune file

Dual filament design reduces maintenance intervals



# High Performance Collision Cell Design

Linear acceleration design is optimized for high speed performance without ion ghosting or cross-talk

MRM speed to 500 MRMs/sec allows determination of more compounds per ion group


High sensitivity with wide mass bandwidth eliminates the need to “tune on your compound” for optimum sensitivity

**“Helium Quenching”** neutral noise reduction increases S/N




# Collision Cell Details

(Cell housing and 2 rods removed for clarity)



# Technical Advantage: Collision Cell



Highly energetic helium metastables are produced in an EI source and can pass through the electrical fields unaffected – producing “neutral noise” at the HED-EM.


## **REMEMBER:**

In GC/MS, Chemical Noise is much higher than Neutral Noise.

In GC/MS/MS with the elimination of Chemical Noise, Neutral Noise is a critical source of noise.



# Technical Advantage: Collision Cell Quench Gas








# Effect of “Helium Quench”

tested with  
hexachlorobenzene

## Noise Comparison





# High Performance Ramped-Iris Detector


Ultra low neutrals noise

Long life and high linearity

“Gain Normalization” corrects tune file for detector aging to allow repeatable long term method sensitivity



# Ramped-Iris Detector



# Agilent 7000B GC/MS/MS

Outstanding sensitivity 100fg of OFN on column at 500:1 S/N RMS in MS/MS mode using AUTOTUNE parameters verified at customer installation.

1050 amu Mass Range

500 MRM/sec Speed

Reliable heated monolithic gold plated hyperbolic quadrupoles

Differentially pumped vacuum system

Helium quenching collision cell technology

Agilent 7890 GC with Capillary Flow technology.

MassHunter Software



# 7000B EI MRM Specifications

## EI MRM Spec:

- S/N 500:1 on 100fg OFN of the 272:222 or 272:241 transition
  - Installation test: as 100fg; splitless
- 8% RSD on peak area of the 272:222 or 272:241 transition



# 7000B CI Specifications

## PCI MRM Spec:

- S/N 50:1 on 100fg BZP MRM of 183:105 transition

## NCI SIM Spec:

- S/N 200:1 on 10fg OFN MS1 SIM of m/z 272



# Autotune

Proprietary program tunes the source, mass analyzer, and detector for (as applicable):

- Ion transmission
- Mass axis calibration
- Mass resolution
- Detector gain vs Voltage

Autotune settings are saved with the method for repeatable method performance.

Manual Tune override is available



# MassHunter MS Workstation Software

Modern software interpretation of the proven industry standard GC/MS ChemStation platform

Single software platform for all Agilent MS Systems

- LC/SQ , LC/MS/MS, LC/TOF, LC/QTOF
- GC/SQ, GC/MS/MS, GC/QTOF (future)
- ICP/MS




# New Agilent Multimode Inlet Improvements

Turn-top easy liner exchange

Standard 11mm septa

Air plus CO<sub>2</sub>, N<sub>2</sub>  
cryogenic cooling

No leaks at liner



Standard liner  
dimensions

Standard column nut



# Agilent Multimode Inlet Features

## Hardware

Temperature range of -160 to 450°C

Heating ramps as fast as 900°C/min

Septum/Liner Easily Exchangeable

Injection Modes: Hot S/SL, Cold S/SL, all in pulsed mode, solvent vent mode, residue removal mode

Support for single stroke injections from 0.1 µL to 250 µL

## Software


Ten temperature ramps

Solution for solvent vent timing

Fully integrated into ChemStation, MSD ChemStation, EZChrom, MassHunter




# Large Volume (LVI) Injection of Triazine Herbicides



**1µL Cold Splitless inj**

**1µL Hot Splitless inj**




**1, 2, 3 and 4 µL Cold splitless injections (ethyl acetate)**

# Backflush: Many Advantages for GC/MS Analysis of Complex Samples ('Dirty Matrices')

- Provides more consistent GC retention times
- Provides better, more consistent MS spectra through a sample sequence
  - Reduces **chemical noise** that may increase during a sequence of samples due to small carryover of matrix from sample to sample
  - Higher quality quantitation (no increase in interfering ions during analysis sequence)
- Reduces contamination (and cleaning frequency) for the source
- Reduces analysis time (more samples per day)
- Increases lifetime of analytical column
- Multiple ways to setup the backflush using Capillary Flow Technologies
  - QuickSwap, Deans Switch, 2-Way Splitter with Makeup, 3-Way Splitter with Makeup, and Purged Ultimate Union




# Without Backflush: A Serious Problem




Overlay of two chromatograms of a blank extract injected BEFORE (A) and AFTER (B) three injections without backflush

# Without Backflush: Changes in Retention Time



Overlay of two chromatograms of a blank extract injected BEFORE (A) and AFTER (B) three injections without backflush


# Without Backflush: Increased Background




Overlay of two chromatograms of a blank extract injected BEFORE (A) and AFTER (B) three injections without backflush

# Backflush

## During GC Run



## After GC Run



# With Backflush: No Increased Background (Less Spectral Noise) and Consistent Retention Times

Abundance

4.6e+07  
4.4e+07  
4.2e+07  
4e+07  
3.8e+07  
3.6e+07  
3.4e+07  
3.2e+07  
3e+07  
2.8e+07  
2.6e+07  
2.4e+07  
2.2e+07  
2e+07  
1.8e+07  
1.6e+07  
1.4e+07  
1.2e+07  
1e+07  
8000000  
6000000  
4000000  
2000000  
0

TIC: lettuce\_10\_ppb.D\data.ms


TIC: lettuce\_100\_ppb.D\data.ms

TIC: lettuce\_5\_ppb.D\data.ms

Data provided by user  
in Almeria, Spain

Stable retention times  
and baseline . . . less  
chemical noise

Time



Overlay of three chromatograms of lettuce extract run with 2 min of back flush





# Purged Ultimate Union for GC/MS

## *New, simple CFT technology for backflush*

Vent-less injector and column maintenance (first column)

Improved backflush


- 1/4 the time and/or lower pressure
- Compatible with diffusion pump MSD

Less loss of sensitivity (less flow added = less dilution)

Option of different phase/dimension for each column

Less complex

- Simple GC configuration (existing supplies)
- Setpoints calculated by GC
- Pressure pulsed injections OK
- Constant flow mode easy!!



# Analytical Reality of GC/MS/MS Methods

**GC/MS/MS needs backflush as much or more than GC/MS to avoid “invisible” problems:**

**source contamination  
loss of sensitivity**

**And to avoid the visible problems:  
changing retention times**



# Method Building and Optimization Tools for the 7000 Series Triple Quad GC/MS/MS

1. **Dynamic MRM for GC/MS –  
optimization of MRM method for  
best analytical performance**
2. **Design Experiment Assistant –  
Create MRM parameter  
optimization experiments**
3. **Analyze Experiments Assistant –  
Automatically analyze MRM  
optimization experiments**



# Dynamic MRM Assistant

Agilent MassHunter DynamicMRM Assistant

Open Samples...

Parameters

Optimize

Save Method As

Exit

Help

Samples

- GR-0902-CH-LDD-5uL-082409-1.D
  - Acquired Transitions
  - Optimized Transitions
  - Optimized Acquisition Method


Optimized Transitions

	Time Segment	Compound Name	Transition	RT	Area	Acquired Dwell	Optimized Dwell
	40	Norfurazon	+MRM: 145-> 145	26.30	840.7	15.0	24.6
	40	Desmethylnorfurazon	+MRM: 145-> 95	26.31	855.1	15.0	40.6
	40	Norfurazon	+MRM: 145-> 95	26.31	855.1	15.0	40.6
	41	Endosulfan sulfate	+MRM: 387-> 253	26.66	324.4	15.0	106.1
	41	Endosulfan sulfate	+MRM: 387-> 217	26.67	188.5	15.0	139.4
	41	DDT, p,p'	+MRM: 237-> 165	26.90	7282.5	15.0	54.0
	41	DDT, p,p'	+MRM: 235-> 165	26.90	11759.2	15.0	34.2
	42	Endosulfan sulfate	+MRM: 387-> 253	26.66	324.4	15.0	106.1
	42	Endosulfan sulfate	+MRM: 387-> 217	26.67	188.5	15.0	139.4
	42	DDT, p,p'	+MRM: 237-> 165	26.90	7282.5	15.0	54.0
	42	DDT, p,p'	+MRM: 235-> 165	26.90	11759.2	15.0	34.2
	43	Propiconazole-I	+MRM: 261-> 175	27.05	93.9	15.0	69.4
	43	Propiconazole-I	+MRM: 259-> 173	27.04	247.2	15.0	42.2
	43	Propiconazole-II (CAS # 60207-90-1)	+MRM: 259-> 173	27.04	247.2	15.0	42.2
	43	Terbuconazole	+MRM: 250-> 163	27.35	1300.0	15.0	56.5
	43	Propiconazole-II (CAS # 60207-90-1)	+MRM: 259-> 191	27.05	138.4	15.0	56.4
	43	Terbuconazole	+MRM: 250-> 125	27.35	4906.1	15.0	28.9
	43	Propiconazole-II (CAS # 60207-90-1)	+MRM: 259-> 69	27.04	237.2	15.0	52.5
	43	Imiprothrin	+MRM: 151-> 123	27.08	290.4	15.0	21.8
	43	Imiprothrin	+MRM: 123-> 95	27.08	7978.7	15.0	20.9
	43	Terbuconazole	+MRM: 125-> 89	27.35	8608.7	15.0	21.8
	43	Imiprothrin	+MRM: 123-> 81	27.08	23357.7	15.0	10.0
	44	Terbuconazole	+MRM: 250-> 163	27.35	1300.0	15.0	56.5
	44	Terbuconazole	+MRM: 250-> 125	27.35	4906.1	15.0	28.9
	44	Terbuconazole	+MRM: 125-> 89	27.35	8608.7	15.0	21.8
	45	Iprodione	+MRM: 314-> 271	28.32	596.0	13.0	62.3
	45	Iprodione	+MRM: 314-> 245.05	28.32	10079.1	13.0	15.1
	45	Chrysene-D12	+MRM: 240-> 236	28.28	665511.5	13.0	10.0
	45	Chrysene-D12	+MRM: 240-> 212	28.28	209467.5	13.0	10.0
	46	Fenprothrin	+MRM: 265-> 210	28.91	3509.9	13.0	14.5
	46	Methoxychlor	+MRM: 227-> 169	28.78	19111.1	13.0	10.0


The Dynamic MRM Assistant optimizes dwell times and time segments for sensitivity.



# GC/MS/MS Software Modules



# Instrument Control UI/Navigation



# MS/MS Acquisition Method Editor Dialog

QQQ Method Editor

Tune file  
atunes.tune.xml

Ion source  
Ion source: ESI  
Source temp.: 230 C  
Electron Energy:  
☐ Tune setting eV  
☐ Fixed eV  
☒ Variable by time segment

☐ Run time 1 min  
Solvent delay 3 min

Time filtering  
☒ Peak with 0.07 min

Detector setting  
☐ Gain ☒ Delta EMV

Time segments

#	Time	Scan Type	Polarity	Electron Energy
1	3.0	MRM	Positive	70
2	4.0	MRM	Positive	70

cycle/s ms/cycle

Timed events

Time	Event
------	-------

☒ Display Timed Events

Acquisition

Scan segments

Compound Name	ISTD?	Precursor Ion	MS1 Resolution	Product Ion	MS2 Resolution	Dwell
Compound1	<input type="checkbox"/>	400	Unit	200	Unit	200
Compound 2	<input type="checkbox"/>	350	Unit	200	Unit	200

Single window acquisition setup

OK Apply Reset Help





# Some GC/MS/MS Applications

- Pesticides in fruits/veges/marine biota
- PCBs
- Nitro-PAHs in air particulates
- Organo-tins
- Melamine in milk and baby formula
- Drug Screens
- PBDE Flame Retardants
- Anabolic Steroids in urine/blood

