

Replacing Mercury Reagents in the Lab

Compiled by INFORM, Inc.

June 13, 2003

There are many resources for those who wish to reduce the use of mercury fixatives and mercury preservatives in the laboratory, including:

- Resources related to determining which products contain mercury.
- Links to Listserv postings by histologists who have switched to mercury-free substitutes for B-5 fixative, Zenker's solution, and PVA parasitology transport vials.
- Publications on this topic.
- Contact information for laboratory personnel willing to share their experience on this topic.

Resources Related to Determining Which Products Contain Mercury

- Thimerosal is a commonly-used mercury preservative, but it can be called by many names. See <http://cerhr.niehs.nih.gov/CERHRchems/Thimerosal.html> for a list of Thimerosal synonyms.
- Dana Farber Cancer Center surveyed their vendors to determine which antibody products contained mercury. See:
 - http://research.dfci.harvard.edu/ehs/mercury_at_dfci.htm, which includes a link to a list of companies that use Thimerosal as preservative in HRP solutions and antibodies
 - http://research.dfci.harvard.edu/ehs/think_before_using_the_sink, which includes a list of companies that do not use Thimerosal or other mercury preservatives
 - For more information on the compilation of these lists, contact Melissa K. McCullough, MS, CHMM, Dana Farber Cancer Center, Melissa_McCullough@dfci.harvard.edu, (617) 632-2619.

Histonet Listserv Archives Postings

The Histonet listserver is a listserver for the Histology profession. Many postings appear in the archives where people discuss switching from mercury laboratory chemicals.

The main page for the Histonet Listserve Archives is at

<http://www.histosearch.com/histonet.html>

Below are links to specific posts providing contact information for and opinions of people who have or are switching from mercury fixatives (B-5 and Zenker's):

- Jude Carpenter of FAHC Histopathology in Vermont posted just last January (2003) that they switched to zinc formalin from Anatech after testing 6

possible substitutes.

<http://www.histosearch.com/histonet/Jan03/RE.B5Substitutes.html>

- Bryan R. Hewlett of Hamilton Regional Laboratory Medicine Program in Ontario said in May 2002 that they were testing B-5 substitutes.
<http://www.histosearch.com/histonet/May02/Re.B5B.html>
- Scott Taft surveyed Histonet members on B-5 use and substitutes in October of 2001. Previous postings identified Scott's email address as ss336@yahoo.com. The results of his survey are posted at <http://www.histosearch.com/histonet/Oct01/ResultsofB5fixativesurvey.html>.
- Teri Johnson of the Stowers Institute for Medical Research in Missouri posted in October of 2001 that they had already switched to zinc formalin.
<http://www.histosearch.com/histonet/Oct01/Re.B-5fixativequestions.html>
- Bob Richmond of Tennessee noted in September of 2000 that overnight fixation in formalin solution can replace fixing in mercury B-5.
<http://www.histosearch.com/histonet/Sep00A/Re.EnteringthemercurydispC.html>
- James Hall of University College London (rmkdhjh@ucl.ac.uk) posted in September of 2000 regarding an alternative Masson technique that does not require mercury.
<http://www.histosearch.com/histonet/Sep00A/Massonalternative.html>
- Sheila Tapper of SMDC Health System in Duluth, Minnesota (STapper@smdc.org) posted in August 1999 that they performed a blind study of various substitutes before deciding on B-Plus.
<http://www.histosearch.com/histonet/Aug99A/RE.B-5SubstitutesA.html>
- Cindy Higgerson of Memorial Hospital in Illinois (chiggerson@memhosp.com) posted in August of 1999 that they were switching from B-5 to a mercury-free alternative. <http://www.histosearch.com/histonet/Aug99A/Re.B-5Substitutes.html>
- Tim Morken of the CDC (tim9@cdc.gov) posted a protocol for switching in October of 1998 at <http://www.histosearch.com/histonet/Oct98/Re.B-5FixativeSubstitutesD.html>

Publications Related to Replacing Mercury Reagents

- A review of the merits, and mechanism and literature of zinc formalin: Dapson R., "Fixation for the 1990's," *Biotechnic & Histochemistry*, 1993, 68(2):75-82.

- Comparison of mercury PVA and alternatives: Jensen, B., et al., "Comparison of Polyvinyl Alcohol Fixative with Three Less Hazardous Fixatives for Detection and Identification of Intestinal Parasites," *Journal of Clinical Microbiology*, April 2000, P. 1592-1598, Vol. 38, No. 4.
- Replacing mercury-containing laboratory chemicals in Duluth: Harvie, J., "Eliminating Mercury Use in Hospital Laboratories: A Step toward Zero Discharge," *Public Health Reports*, July/August 1999, 114(4):353-8. This article can be downloaded at <http://noharm.org/details.cfm?type=document&id=256>. The author's contact information has changed since this paper was published to Jamie Harvie, Institute for a Sustainable Future, 32 E. 1st Street, Suite 206, Duluth MN 55802, 218-525-7806, harvie@isfusa.org.
- "Improved Methods for Staining Fibrin", Chuck Churukian, *Journal of Histotechnology*, in press June 2003.
- "The Manual of the Special Stains Laboratory," 2002 Edition - \$35.50, includes postage. Send check/money order to Department of Pathology, c/o Chuck Churukian; University of Rochester Medical Center; Box 626; Rochester NY 14642; You may contact Mr. Churukian at 585-275-1864 (FAX 585-273-1027, Churukian@aol.com)

Laboratories Willing to Share Their Experience

- Pierre Gonyon at Trinity Health/St. Joseph's in Ann Arbor posted an email to the H2E listserve in January of 2002 detailing how his hospital did a blind study comparing the mercury B-5 to the zinc and the histologists couldn't tell the difference. Protocol/product information for switching from B-5 to B-Plus fixative is available from him: Pierre Gonyon, Safety Department, Saint Joseph Mercy Health System, Ann Arbor, MI 48106, (734)712-3315, gonyonpl@trinityhealth.org.
- Lisa Meyer of Fisher-Titus Medical Center (lmeyer@ftmc.com) switched several years ago.
- Crystal Sands of NorDx Laboratories (207-885-7807, sandsc@mmc.org) has been making these changes They replaced;
 - Mercury PVA after conducting a trial of 4 alternatives. In the clinical lab they replaced PVA with Zinc PVA. The lab conducted trials of SAF, Ecofix, Spin-Con and the Zinc PVA.
 - In Histology, they eliminated Zenkers completely. After two years and many product trials, including B Plus, Z5, B5 Fixative Substitute, and AZF, the pathologists agreed to switch from B5 to AZF.
- Deb Goodman, Housekeeping Supervisor at SSM Health Care (Debbie_Goodman@ssmhc.com) has reported that switching was easy, the

reagents did not cost more, and they saved on hazardous waste costs. They did side-by-side trials and retrained their technicians. The new reagents require fewer steps.